





























































































Reusable Dental Diamond Burs

SHAPE	PART #	SHANK Ø	HEAD L	GRIT
CYLINDER AND END CUTTING				
	111-836G-010 111-836-010	R 010	4 mm	
	111-837G-012	R 012	7 mm	
CYLINDRICAL, ROUND END				
	111-881G-010	R 010	8 mm	
	111-881SG-012	R 012 S	8 mm	
EGG				
	111-833SG-023	R 023	4.2 mm	
	111-833EF-023	R 023	4.2 mm	
FINISHING, TAPERED-POINTED				
	111-859EF-012	R 014	10 mm	
FLAME				
	111-862-012	R 012	6 mm	
	111-862G-012 111-862F-012 111-862EF-012	R 010	8 mm	
	111-863G-012	R 008	10 mm	
FOOTBALL				
	111-831F-016	R 016	3.5 mm	
	111-830G-021	R 021	5 mm	
	111-830SG-023 111-830F-023 111-830EF-023	R 023	5 mm	
INVERTED CONICAL W/ COLLAR				
	111-806G-014	R 014	2.2 mm	
OCCUSAL-GRINDER				
	111-811SG-033	R 033	4.5 mm	
	111-811LSG-016	R 037	7 mm	
PEAR				
	111-809-008	R 008	1.6 mm	
	111-808-012	R 012	2.7 mm	
	111-808G-012	R 012	2.7 mm	
ROUND				
	111-801G-023	R 023	2 mm	
TAPERED, FLAT END				
	111-847SG-014	R 014	8 mm	

SHAPE	PART #	SHANK Ø	HEAD L	GRIT
TAPERED, FLAT END (CONT.)				
	111-847SG-016	R 016	8 mm	
	111-848SG-016	R 016	10 mm	
	111-847SG-018	R 018	8 mm	
TAPERED, POINTED				
	111-890F-008	R 008	3 mm	
	111-890LF-008	R 008	3 mm	
	111-858G-014	R 014	8 mm	
	111-859G-014 111-859F-014	R 014	10 mm	
TAPERED, ROUND EDGE				
	111-847RSG-016	R 016	8 mm	
TAPERED, ROUND END				
	111-191G-012	R 012	7 mm	
	111-852G-012	R 012	10 mm	
	111-194G-012	R 014	8 mm	
	111-852G-014	R 014	10 mm	
	111-191SG-Q-016	Q 016	7 mm	
	111-194SG-016	R 016	8 mm	
	111-852SG-016 111-852-016	R 016	10 mm	
	111-191SG-018	R 018	7 mm	
	111-852SG-018	R 018	10 mm	
	111-194SG-021	R 021	8 mm	
	111-191SG-025	R 025	7 mm	
TORPEDO, CYLINDRICAL				
	111-868G-012	R 012 S	8 mm	
	111-885G-012	R 012	8 mm	
TORPEDO, TAPERED				
	111-878G-014	R 014	8 mm	
	111-878G-016	R 016	8 mm	
	111-878SG-018	R 018	8 mm	
WHEEL				
	111-909G-040	R 040	1.5 mm	

Grit Size:  Super Coarse 150-180 µm  Coarse 125-150 µm  Medium 90-106 µm  Fine 38-45 µm  Extra Fine 20-30 µm

Shank Type: Q = Friction Grip, Short, 16 mm / R = Friction Grip, 19 mm